

KYRA® 500 Stirrup

Strong. More Padding.



The KYRA® Comfort™ 500 stirrup is designed for patients with a BMI of greater than 29 with a large flexible lateral fin and can support patients up to 500 lbs/227 kg weight capacity.

KYRA boots will safely accept larger legs with 1" of supportive boot padding.

All Comfort 500's come standard with the durable Secure-Lok™ straps which help reduce infection control and the patented CoreGrip™ technology which makes the stirrup stronger and lighter vs other stirrup brands.

Specifications

- KYRA® Comfort 500 weight: 14.6 lbs/6.6 kg each
- Lightweight and strong with KYRA's CoreGrip™ technology
- More working space with wider abduction range and boot rotational angle
- Lithotomy Range: +85° to -55°, Abduction Range: +25° to -9°, Boot Rotational Angle: +/-25°

KYRA5100 KYRA 500 Stirrups w/ Secure-Lok™ straps

Purple is the New Yellow



ADVOCATES
FOR PATIENT
SAFETY AND
QUALITY
EQUIPMENT