



QUIK VAC

BEAN BAG POSITIONERS



QV11
400mm x 350mm



QV30
730mm x 900mm



QV31
980mm x 900mm



QV32
1200mm x 900mm



QV33
1360mm x 900mm

Quik Vac Bean Bag Positioners are a U shaped positioner and available in Head and Neck plus optional widths to accommodate varying procedures and sizes.

Features

- Suction used to stabilise and retain bean bag shape that holds the patient securely
- By opening the valve, the bean bag softens and patient can be repositioned
- Safe method for positioning and stabilising patients undergoing procedures
- Various sizes to accommodate different procedures and weights
- Easily wiped over with PH Neutral detergent or wipe
- They conform comfortably to the patients shape
- 1 Year manufacturers warranty
- Latex FREE

Code Dimensions

QV11	400mm (w) x 350mm (l)
QV30	730mm (w) x 900mm (l)
QV31	980mm (w) x 900mm (l)
QV32	1200mm(w) x 900mm (l)
QV33	1360mm(w) x 900mm (l)

Our focus is patient care



ADVOCATES
FOR PATIENT
SAFETY AND
QUALITY
EQUIPMENT