

KYRA[®] Armboards




KYRA Armboard, Swivel, Large KYRA2300

Allows for effective positioning of the arms to help prevent potential patient injury. The hand rest is designed to enable the surgeon to place the patient's hand in the desired orientation. Pad included.

Armboard Weight: 2.8 kg
Armboard Post Size: 46.2 cm x 31.7 cm x 1.6 cm




Included



KYRA2340
Armboard, Swivel
x1 Included

Required



KYRA4800
Round Post Rail Clamp
x1 Required


KYRA Armboard, Swivel, Small KYRA2304

Allows for effective positioning of the arms to help prevent potential patient injury. The hand rest is designed to enable the surgeon to place the patient's hand in the desired orientation. Pad included.

Armboard Weight: 2.7 kg
Armboard Post Size: 40.6 cm x 16.5 cm x 1.6 cm




Included



KYRA2340
Armboard, Swivel
x1 Included

Required



KYRA4800
Round Post Rail Clamp
x1 Required

The Original. Without Compromise.



ADVOCATES
FOR PATIENT
SAFETY AND
QUALITY
EQUIPMENT

KYRA® Armboards and Pads



KYRA® Trigger Armboard

Integrated clamp snap locks onto rail at any position. Designed to prevent fluid ingress. Radiolucent, ultra durable strong material. Eliminates gear grinding and pinch hazards. One-hand adjustable. Pad not included

Dimensions: 66 cm x 14 cm
Armboard Weight: 1.8 kg
Patient Weight Capacity: 227 kg

Specify Rail Type

KYRA2000 US
KYRA2005 EU
KYRA2006 UK
KYRA2007 DY
KYRA2008 JP



KYRA® Comfort™ Armboard Pad

Traditional pads with soft pressure relieving foam. Available in various thicknesses. Fits armboards 66 cm x 14 cm. Includes 2 Soft Armboard Straps.

Standard Pad Weight:
KYRA2009/2010 0.36 kg
KYRA2015/2016 0.41 kg

Specify Pad Size

KYRA2009 2.5 cm
KYRA2010 5 cm
KYRA2015 7.6 cm
KYRA2016 10 cm



KYRA® Clamshell Armboard Pad

Integrated clamshell strap safely encapsulates arm. Available in various thicknesses. Fits armboards 66 cm x 14 cm.

Clamshell Padded Area: 25.4 cm

Clamshell Padded Weight:
KYRA2050/2051 0.43 kg
KYRA2052/2053 0.45 kg

Specify Pad Size

KYRA2050 2.5 cm
KYRA2051 5 cm
YRA2052 7.6 cm
KYRA2053 10 cm



The Original. Without Compromise.



ADVOCATES
FOR PATIENT
SAFETY AND
QUALITY
EQUIPMENT