

KYRA Curve™ Spine Frame System

KYRA9000K



Provides a stable platform for spine surgeries including laminectomies, decompressions and disc surgery. Intraoperatively adjustable to the desired height to get optimal surgical site exposure. C-arm compatible. Fits spine table rails and tabletops. Unique silicone straps provides better cleaning vs. hook & loop straps. Pad set included.

Patient Weight Capacity: 272 kg
Device Weight: 14.24 kg



- Integrated gel pads
- Allows for more abdomen space



- Durable silicone straps
- Hook & loop free

Included



KYRA9000
Curve Spine Frame
x1 Included

KYRA9010
Curve Crank Handle
x1 Included

KYRA1100
Curve Spine Frame Cart
x1 Included

Optional



KYRA9020
KYRA Curve Comfort Cover Set
x6 pack

The Original. Without Compromise.



ADVOCATES
FOR PATIENT
SAFETY AND
QUALITY
EQUIPMENT

KYRA Curve™ Spine Frame System Accessories



KYRA Curve Spine Frame Replacement Pads

KYRA6450-P

Pair of replacement pads for the Curve Spine Frame.
No special tools necessary to remove or attach the pads.



KYRA Curve Spine Frame Cart KYRA1100

Easy to push/pull with ergonomic handles. Cart keeps KYRA Curve Spine Frame secured during transport to prevent damage. Locking casters.

Cart Height: 103.5 cm
Base Width: 54.93 cm x 54.93 cm
Cart Weight: 16.1 kg



KYRA Curve Crank Handle KYRA9010

Rotational handle for the Curve Spine Frame to create lordosis. Laser marked and windowed to help ensure positive connection with the Curve Spine Frame.



The Original. Without Compromise.



ADVOCATES
FOR PATIENT
SAFETY AND
QUALITY
EQUIPMENT