


All KYRA® Stirrups provide unparalleled pelvic site access, larger boot sizing and better infection control.



Curved rod provides more ROM and working space

- The unique patent-pending design provides 22° more lithotomy Range of motion than traditional stirrups
- Great for robotic surgery
- Design protects the gas spring for longer life



Larger boot and more padding for today's patients

- Fits patient's legs into the boots, not on top of the boot
- No exposed hard plastic edges
- Flexible lateral fin



VELCRO® free straps

- Medical grade silicone straps are easy to clean and snap-in securely
- Strong but flexible
- Non-propagating material



The Original. Without Compromise.



ADVOCATES
FOR PATIENT
SAFETY AND
QUALITY
EQUIPMENT

The KYRA® Comfort™ Family of Stirrups

- Designed to improve infection control with optional Secure-Lok™ straps
- Lithotomy Range: +85° to -55°, Abduction Range: +25° to -9°, Boot Rotational Angle: +/-25°
- Boots sized to fit today's larger patient anatomy
- Easy to clean durable silicone straps
- More pressure reducing padding
- Lightweight and strong



Included

OR

KYRA3000
Comfort 350
Stirrups w/ hook
& loop straps
x1 Pair Included

KYRA3100
Comfort 350
Stirrups w/ Secure
Lok™ silicone straps
x1 Pair Included

Required

KYRA4100
Secure-Release™ Clamp
x2 Required

Optional

KYRA1100
Curve Spine
Frame Cart



Included

OR

KYRA5000
Comfort 500
Stirrups w/ hook
& loop straps
x1 Pair Included

KYRA5100
Comfort 500
Stirrups w/ Secure
Lok™ silicone straps
x1 Pair Included

Required

KYRA4100
Secure-Release™ Clamp
x2 Required

Optional

KYRA1100
Curve Spine
Frame Cart



Included

KYRA8100
Comfort 800 Stirrups
w/ Secure-Lok™
silicone straps
x1 Pair Included

Required

KYRA4100
Secure-Release™ Clamp
x2 Required

Optional

KYRA1100
Curve Spine
Frame Cart



KYRA® Single™ Stirrup Boot Covers

Infection Control for your stirrups

- Protect against cross contamination
- Fits Yellow and Purple branded stirrups
- Single use, protective barrier
- For infectious patients (known or suspected)
- For patients with open wounds
- Decrease room turnover time

KYRA8149 KYRA Single Boot Covers for KYRA 350/500/800 Stirrups (10 pack)
KYRA8198 KYRA Single Boot Covers for Yellofins® Stirrups (10 pack)

The Original. Without Compromise.



ADVOCATES
FOR PATIENT
SAFETY AND
QUALITY
EQUIPMENT

KYRA® Stirrup Replacement Pads

KYRA® Quality for non-Kyra stirrups.



KYRA6400-P
Fits Yellow stirrups



KYRA6420-P
Fits Blue non-fin & Basic stirrups



KYRA6410-P
Fits Bariatric stirrups



- Kyra® Medical® offers true replacement pads for yellow stirrups, blue non-fin stirrups and bariatric stirrups
- Sold in pairs, our replacement pads match the quality of the original equipment
- Call your local KYRA® sales representative today and book a stirrup pad assessment.
Your rep will provide you with a no-cost evaluation of your stirrup's pads
- If your stirrups need more than just new pads, ask for a demo of Kyra® Stirrups – The Purple Stirrups!

KYRA6400-P Stirrup pads, fits Yellow stirrups

KYRA6420-P Stirrup pads, fits Blue non-fin and Basic stirrups

KYRA6410-P Stirrup pads, fits Bariatric stirrups

The Original. Without Compromise.

1300 655 160 | Info@Endovision.com.au | Endovision.com.au

Endovision

ADVOCATES
FOR PATIENT
SAFETY AND
QUALITY
EQUIPMENT