

Evolution Patient Positioners



VascoCare Medical has introduced innovative new Supine Positioners to the EVOLUTION Patient Positioner Range. These provide enhanced pressure relief with a proven anatomical design.




Supine Head Positioner

The head positioner has a unique design that redistributes head weight and provides better pressure relief than regular head rings. It is made of two components - comfortable visco-elastic foam on the inside and a strong, seamless outer made of a durable hi-tech membrane. Suitable for ENT and ophthalmic surgery.




Supine Knee Positioner

The leg and knee positioner has a unique design that supports the leg and knee, offering enhanced pressure relief compared to regular positioners. It is made of two components - comfortable visco-elastic foam on the inside and a strong, seamless outer made of a durable hi-tech membrane.



Supine Foot Positioner

The foot positioner has a unique design that redistributes weight and provides enhanced pressure relief on the calcaneus bone compared to regular heel pads. It is made of two components - comfortable visco-elastic foam on the inside and a strong, seamless outer made of a durable hi-tech membrane.



Key Features

- Unique and innovative design for enhanced pressure relief and weight redistribution
- Made of two components: comfortable Visco-elastic foam on the inside and a strong, seamless Dual Coated Hi-Tech Membrane on the outside
- Outer coating is highly durable, with excellent abrasion resistance and puncture/burst resistance
- Flexible design improves conformity to the skin
- Provides zero interface pressure under bony prominences
- Strong bond between foam and coating ensures a long-lasting product
- Outer coating acts as a fluid-proof barrier, preventing bacterial growth
- Compatible with CT Scanners without causing artifacts.
- Head positioners offer superior neck support and head stability compared to regular head rings.
- Knee positioners provide excellent leg support compared to regular domed positioners.
- Foot positioners offer exceptional foot support and zero pressure under the heel, unlike regular heel pads.
- Competitively priced to be a cost-effective alternative to regular positioners.
- Latex-free.
- Comes with a two-year product warranty.

VSHP23 Supine Head Positioner Medium Style 2

VSHP24 Supine Head Positioner Large Style 2

VSKP13 Supine Knee Positioner Medium

VSFP13 Supine Foot Positioner Medium

Our focus is staff and patient safety



ADVOCATES
FOR PATIENT
SAFETY AND
QUALITY
EQUIPMENT

Quik Vac Bean Bag Positioners



Quik Vac Bean Bag Positioners are a U shaped positioner and available in Head and Neck plus optional widths to accommodate varying procedures and sizes.



QV11
400mm x 350mm



QV30
730mm x 900mm



QV31
980mm x 900mm



QV32
1200mm x 900mm



QV33
1360mm x 900mm

Key Features

- Suction used to stabilise and retain bean bag shape that holds the patient securely
- By opening the valve, the bean bag softens and patient can be repositioned
- Safe method for positioning and stabilising patients undergoing procedures
- Various sizes to accommodate different procedures and weights
- Easily wiped over with PH Neutral detergent or wipe
- They conform comfortably to the patients shape
- 1 Year manufacturers warranty
- Latex FREE.

- QV11** 400mm (w) x 350mm (l)
- QV30** 730mm (w) x 900mm (l)
- QV31** 980mm (w) x 900mm (l)
- QV32** 1200mm (w) x 900mm (l)
- QV33** 1360mm (w) x 900mm (l)

Our focus is staff and patient safety



ADVOCATES
FOR PATIENT
SAFETY AND
QUALITY
EQUIPMENT