



*Purple is the New Yellow.*



ADVOCATES  
FOR PATIENT  
SAFETY AND  
QUALITY  
EQUIPMENT

[Endovision.com.au](http://Endovision.com.au)

# Reduce your downtime with KYRA<sup>®</sup> Medical<sup>®</sup> Onsite Repair Kits!

Cost effective stirrup maintenance for Today's Modern OR<sup>™</sup>!



## No Repair Cost Surprises

KYRA<sup>®</sup> Medical<sup>®</sup> publishes price lists for ALL anticipated stirrup repairs<sup>1</sup>. You know EXACTLY what you will pay for any repair.



## Low-Cost No-Hassle Training

KYRA<sup>®</sup> Medical<sup>®</sup> provides online training that will certify your biomedical/facility teams at no cost to your facility.



## Onsite Repair Kits Reduce Maintenance Costs

KYRA<sup>®</sup> Medical<sup>®</sup> provides reasonably priced repair kits that allows you to make onsite stirrup repairs 90%<sup>2</sup> of the time.



## Reduce Costly Equipment Downtime

Only KYRA<sup>®</sup> Medical<sup>®</sup> allows your biomedical team to make most repairs on site; reducing costly equipment down time.



Please contact us for Repair Kit pricing

1. In addition, Customer pays inbound and outbound shipping charges; cost of repairs are equal to repair kits X 20% for labor. 2. These repairs require custom jigs and tools and are not available as kits.

*The Original. Without Compromise.*

**Endovision**

ADVOCATES  
FOR PATIENT  
SAFETY AND  
QUALITY  
EQUIPMENT

1300 655 160 | [Info@Endovision.com.au](mailto:Info@Endovision.com.au) | [Endovision.com.au](http://Endovision.com.au)

All KYRA® Stirrups provide unparalleled pelvic site access, larger boot sizing and better infection control.



## Curved rod provides more ROM and working space

- The unique patent-pending design provides 22° more lithotomy Range of motion than traditional stirrups
- Great for robotic surgery
- Design protects the gas spring for longer life



## Larger boot and more padding for today's patients

- Fits patient's legs into the boots, not on top of the boot
- No exposed hard plastic edges
- Flexible lateral fin



## VELCRO® free straps

- Medical grade silicone straps are easy to clean and snap-in securely
- Strong but flexible
- Non-propagating material



*The Original. Without Compromise.*



ADVOCATES  
FOR PATIENT  
SAFETY AND  
QUALITY  
EQUIPMENT

# The KYRA® Comfort™ Family of Stirrups


- Designed to improve infection control with optional Secure-Lok™ straps
- Lithotomy Range: +85° to -55°, Abduction Range: +25° to -9°, Boot Rotational Angle: +/-25°
- Boots sized to fit today's larger patient anatomy
- Easy to clean durable silicone straps
- More pressure reducing padding
- Lightweight and strong




**KYRA®  
Comfort™  
350 Stirrups**

Stirrup Weight 13.2 lbs/6 kg  
Patient Weight Capacity 350 lbs/159 kg

**Included**



**OR**

**KYRA3000**  
Comfort 350  
Stirrups w/ hook  
& loop straps  
x1 Pair Included


**KYRA3100**  
Comfort 350  
Stirrups w/ Secure  
Lok™ silicone straps  
x1 Pair Included

**Required**



**Secure-Release™ Clamp**  
x2 Required

**Optional**



**KYRA1000**  
The Purple  
Stirrup Cart



**KYRA®  
Comfort™  
500 Stirrups**

Stirrup Weight 14.6 lbs/6.6 kg  
Patient Weight Capacity 500 lbs/227 kg

**Included**



**OR**

**KYRA5000**  
Comfort 500  
Stirrups w/ hook  
& loop straps  
x1 Pair Included


**KYRA5100**  
Comfort 500  
Stirrups w/ Secure  
Lok™ silicone straps  
x1 Pair Included

**Required**



**Secure-Release™ Clamp**  
x2 Required

**Optional**




**KYRA1000**  
The Purple  
Stirrup Cart



**KYRA®  
Comfort™  
800 Stirrups**


Stirrup Weight 14.6 lbs/6.6 kg  
Patient Weight Capacity 800 lbs/363 kg

**Included**



**KYRA8100**  
Comfort 800 Stirrups  
w/ Secure-Lok™  
silicone straps  
x1 Pair Included

**Required**



**Secure-Release™ Clamp**  
x2 Required

**Optional**



**KYRA1000**  
The Purple  
Stirrup Cart



## KYRA® Single™ Stirrup Boot Covers

*Infection Control for your stirrups*

- Protect against cross contamination
- Fits Yellow and Purple branded stirrups
- Single use, protective barrier
- For infectious patients (known or suspected)
- For patients with open wounds
- Decrease room turnover time

KYRA8149 KYRA Single Boot Covers for KYRA 350/500/800 Stirrups (10 pack)  
KYRA8198 KYRA Single Boot Covers for Yellofins® Stirrups (10 pack)

The Original. Without Compromise.



1300 655 160 | Info@Endovision.com.au | Endovision.com.au

ADVOCATES  
FOR PATIENT  
SAFETY AND  
QUALITY  
EQUIPMENT

# KYRA® Stirrup Replacement Pads

KYRA® Quality for non-Kyra stirrups.



**KYRA6400-P**  
*Fits Yellow stirrups*



**KYRA6420-P**  
*Fits Blue non-fin & Basic stirrups*



**KYRA6410-P**  
*Fits Bariatric stirrups*



- Kyra® Medical® offers true replacement pads for yellow stirrups, blue non-fin stirrups and bariatric stirrups
- Sold in pairs, our replacement pads match the quality of the original equipment
- Call your local KYRA® sales representative today and book a stirrup pad assessment.  
Your rep will provide you with a no-cost evaluation of your stirrup's pads
- If your stirrups need more than just new pads, ask for a demo of Kyra® Stirrups – The Purple Stirrups!

KYRA6400-P Stirrup pads, fits Yellow stirrups

KYRA6420-P Stirrup pads, fits Blue non-fin and Basic stirrups

KYRA6410-P Stirrup pads, fits Bariatric stirrups

*The Original. Without Compromise.*

**Endovision**

ADVOCATES  
FOR PATIENT  
SAFETY AND  
QUALITY  
EQUIPMENT

1300 655 160 | [Info@Endovision.com.au](mailto:Info@Endovision.com.au) | [Endovision.com.au](http://Endovision.com.au)

# KYRA® Armboards




## KYRA Armboard, Swivel, Large KYRA2300

Allows for effective positioning of the arms to help prevent potential patient injury. The hand rest is designed to enable the surgeon to place the patient's hand in the desired orientation. Pad included.

Armboard Weight: 2.8 kg


Armboard Post Size: 46.2 cm x 31.7 cm x 1.6 cm

Included



KYRA2300  
Armboard, Swivel  
x1 Included

Required



Round Post Rail Clamp  
x1 Required




## KYRA Armboard, Swivel, Small KYRA2304

Allows for effective positioning of the arms to help prevent potential patient injury. The hand rest is designed to enable the surgeon to place the patient's hand in the desired orientation. Pad included.

Armboard Weight: 2.7 kg


Armboard Post Size: 40.6 cm x 16.5 cm x 1.6 cm

Included



KYRA2304  
Armboard, Swivel  
x1 Included

Required



Round Post Rail Clamp  
x1 Required



*The Original. Without Compromise.*



ADVOCATES  
FOR PATIENT  
SAFETY AND  
QUALITY  
EQUIPMENT

# KYRA® Armboards and Pads



## KYRA® Trigger Armboard

Integrated clamp snap locks onto rail at any position. Designed to prevent fluid ingress. Radiolucent, ultra durable strong material. Eliminates gear grinding and pinch hazards. One-hand adjustable. Pad not included

Dimensions: 66 cm x 14 cm  
Armboard Weight: 1.8 kg  
Patient Weight Capacity: 227 kg

### Specify Rail Type

KYRA2000 US  
KYRA2005 EU  
KYRA2006 UK  
KYRA2007 DY  
KYRA2008 JP



## KYRA® Comfort™ Armboard Pad

Traditional pads with soft pressure relieving foam. Available in various thicknesses. Fits armboards 66 cm x 14 cm. Includes 2 Soft Armboard Straps.

Standard Pad Weight:  
KYRA2009/2010 0.36 kg  
KYRA2015/2016 0.41 kg

### Specify Pad Size

KYRA2009 2.5 cm  
KYRA2010 5 cm  
KYRA2015 7.6 cm  
KYRA2016 10 cm



## KYRA® Clamshell Armboard Pad

Integrated clamshell strap safely encapsulates arm. Available in various thicknesses. Fits armboards 66 cm x 14 cm.

Clamshell Padded Area: 25.4 cm

Clamshell Padded Weight:  
KYRA2050/2051 0.43 kg  
KYRA2052/2053 0.45 kg

### Specify Pad Size

KYRA2050 2.5 cm  
KYRA2051 5 cm  
YRA2052 7.6 cm  
KYRA2053 10 cm



*The Original. Without Compromise.*



ADVOCATES  
FOR PATIENT  
SAFETY AND  
QUALITY  
EQUIPMENT

# KYRA® Clamps



## Secure-Release™ Clamp

Mounts and locks anywhere on standard OR table rail.  
Prevents false rail engagement. Push button safely releases clamp. Autoclavable.

### Specify Rail Type

KYRA4100 US  
KYRA4200 EU  
KYRA4300 UK  
KYRA4400 DY  
KYRA4500 JP



## Simple Blade Clamp

Mounts and locks anywhere on standard OR table rail ends or at rail notches. A cost effective solution.

### Specify Rail Type

KYRA4700 US  
KYRA4710 EU/UK  
KYRA4730 DY  
KYRA4740 JP



## Accessory Clamp - Round & Rectangular

Mounts and locks anywhere on standard OR table rail.  
Accepts round post up to 1.8 cm and rectangular post 2.5 cm x 0.6 cm. Autoclavable.

### Specify Rail Type

KYRA4800 US  
KYRA4810 EU  
KYRA4820 UK  
KYRA4830 DY  
KYRA4840 JP



*The Original. Without Compromise.*



1300 655 160 | [Info@Endovision.com.au](mailto:Info@Endovision.com.au) | [Endovision.com.au](http://Endovision.com.au)

ADVOCATES  
FOR PATIENT  
SAFETY AND  
QUALITY  
EQUIPMENT

# KYRA Curve™ Spine Frame System

KYRA9000K



Provides a stable platform for spine surgeries including laminectomies, decompressions and disc surgery. Intraoperatively adjustable to the desired height to get optimal surgical site exposure. C-arm compatible. Fits spine table rails and tabletops. Unique silicone straps provides better cleaning vs. hook & loop straps. Pad set included.

Patient Weight Capacity: 272 kg  
Device Weight: 14.24 kg



- Integrated gel pads
- Allows for more abdomen space



- Durable silicone straps
- Hook & loop free

## Included



**KYRA9000**  
Curve Spine Frame  
x1 Included

**KYRA9010**  
Curve Crank Handle  
x1 Included

**KYRA1100**  
Curve Spine Frame Cart  
x1 Included

## Optional



**KYRA9020**  
KYRA Curve Comfort Cover Set  
x6 pack

*The Original. Without Compromise.*



ADVOCATES  
FOR PATIENT  
SAFETY AND  
QUALITY  
EQUIPMENT

# KYRA Curve™ Spine Frame System Accessories



## KYRA Curve Spine Frame Replacement Pads

### KYRA6450-P

Pair of replacement pads for the Curve Spine Frame.  
No special tools necessary to remove or attach the pads.



## KYRA Curve Spine Frame Cart KYRA1100

Easy to push/pull with ergonomic handles. Cart keeps KYRA Curve Spine Frame secured during transport to prevent damage. Locking casters.

Cart Height: 103.5 cm  
Base Width: 54.93 cm x 54.93 cm  
Cart Weight: 16.1 kg



## KYRA Curve Crank Handle KYRA9010

Rotational handle for the Curve Spine Frame to create lordosis. Laser marked and windowed to help ensure positive connection with the Curve Spine Frame.



*The Original. Without Compromise.*

**Endovision**

ADVOCATES  
FOR PATIENT  
SAFETY AND  
QUALITY  
EQUIPMENT

1300 655 160 | [Info@Endovision.com.au](mailto:Info@Endovision.com.au) | [Endovision.com.au](http://Endovision.com.au)



## Lateral Positioning Set KYRA9100K

Secure and comfortable positioning for procedures requiring total lateral stability. Fully adjustable to any angle. Swivel heads conform to anatomical contours of patient and can be securely locked for precise positioning. No disposable needed. The 3-piece set includes 2 Rectangular and 1 Circular Positioners.

Patient Weight Capacity: 272 kg



### Included



**KYRA9120**  
Rectangular Lateral Positioner  
x2 Included



**KYRA9110**  
Circular Lateral Positioner  
x1 Included

### Required



**Round Post Rail Clamp**  
x1 Required

## Secondary Shoulder Supports KYRA2340

Designed to provide non-contact support for extreme Trendelenburg positioning. Flexible swivel adjustment for optimal positioning. Pads included. Sold in pairs.

Patient Weight Capacity: 272 kg

### Included



**KYRA2340**  
Shoulder Supports  
x1 Pair Included

### Required



**Round Post Rail Clamp**  
x1 Required



*The Original. Without Compromise.*



ADVOCATES  
FOR PATIENT  
SAFETY AND  
QUALITY  
EQUIPMENT



*Purple is the New Yellow.*

Endovision

ADVOCATES  
FOR PATIENT  
SAFETY AND  
QUALITY  
EQUIPMENT

[Endovision.com.au](http://Endovision.com.au)