

DISPOSABLE BLADED TROCAR



Disposable laparoscopic trocar kit, the safety you can feel!

CE 2292

Product features



Appropriate Sealing offers good air tightness.



Threaded Design control insertion speed and depth.



Sharp linear blade allows low insertion force and reduced fascial defect.



Smooth conical aiming facilitates easier insertion.



Simple seal detachment for versatility during specimen removal.



Safety Lock indicates when the blade is shielded and exposed.

Specifications

No.	Ref No.	Description
1	VC-100-1-5	bladed, 5*100mm
2	VC-100-1-10	bladed, 10*100mm
3	VC-100-1-12	bladed, 12*100mm
4	VCR-100-2-A	bladed, 2X5mm with smoke filter
5	VCR-100-2-B	bladed, 2X10mm with smoke filter
6	VCR-100-2-C	bladed, 2X12mm with smoke filter
7	VCR-100-3-A	bladed, 1X5mm,2x10mm with smoke filter
8	VCR-100-4-A	bladed, 2X5mm,2x10mm with smoke filter
9	VCR-100-4-B	bladed, 2X5mm,1x10mm, 1x12mm with smoke filter
10	VCR-100-4-C	bladed, 2X5mm,2x12mm with smoke filter

Note: For more product information, please consult our marketing department.

DISPOSABLE BLADELESS TROCAR



Disposable laparoscopic trocar kit, the safety you can feel!

CE 2292

■ Product features



Universal Sealing increases trocar insertion resistance without using reducer.



Large aiming tapered surface facilitates one-handed operation and easy insert working channel.




Threaded Design enables to control puncture speed and depth and fixed firmly that avoids damage to abdomen.



Translucent Design for easier observation of instrument insertion/removal and extraction of specimen.



Pyramid Tip protects inter organs when puncturing.



■ Specifications

No.	Ref No.	Description
1	IIC-100-1-5	Bladeless, 5*100mm
2	IIC-100-1-10	Bladeless, 10*100mm
3	IIC-100-1-12	Bladeless, 12*100mm
4	IIC-100-1-15	Bladeless, 15*100mm
5	IIC-150-1-5	Bladeless, 5*150mm, For Bariatric Surgery
6	IIC-150-1-10	Bladeless, 10*150mm, For Bariatric Surgery
7	IIC-150-1-12	Bladeless, 12*150mm, For Bariatric Surgery
8	IICR-100-2-A	Bladeless 2x5mm with smoke filter
9	IICR-100-2-B	Bladeless 2x10mm with smoke filter
10	IICR-100-2-C	Bladeless 2x12mm with smoke filter
11	IICR-100-3-A	Bladeless 1x5mm, 2x10mm with smoke filter
12	IICR-100-4-A	Bladeless 2x5mm, 2x10mm with smoke filter
13	IICR-100-4-B	Bladeless 2x5mm, 1x10mm, 1x12mm with smoke filter
14	IICR-100-4-C	Bladeless 2x5mm, 2x12mm with smoke filte

DISPOSABLE OPTICAL TROCAR



Disposable laparoscopic trocar kit, the safety you can feel!

CE 2292

■ Product features




Appropriate Sealing
offers good air tightness.



Threaded Design control
insertion speed and depth.



Visualization tip allows
excellent identification.



Smooth conical aiming
facilitates easier insertion.



Simple seal detachment
for versatility during specimen
removal.



Fixed button with advanced
fixation keep scope from moving
during procedures.

■ Specifications

No.	Ref No.	Description
1	IIV-100-1-5	Optical, 5*100mm
2	IIV-100-1-10	Optical, 10*100mm
3	IIV-100-1-12	Optical, 12*100mm
4	IIVR-100-2-A	Optical 2x5mm with smoke filter
5	IIVR-100-2-B	Optical 2x10mm with smoke filter
6	IIVR-100-2-C	Optical 2x12mm with smoke filter
7	IIVR-100-3-A	Optical 1x5mm, 2x10mm with smoke filter
8	IIVR-100-4-A	Optical 2x5mm, 2x10mm with smoke filter
9	IIVR-100-4-B	Optical 2x5mm, 1x10mm, 1x12mm with smoke filter
10	IIVR-100-4-C	Optical 2x5mm, 2x12mm with smoke filter